


# **Lipid emulsion with Propofol**

## **Direct and accelerated emulsion stability and droplet size distribution**

*2020*

LUM GmbH, Berlin, Germany



# Propofol



- Propofol, is a short-acting medication that results in a decreased level of consciousness and lack of memory for events. Its uses include the starting and maintenance of general anesthesia, sedation for mechanically ventilated adults, and procedural sedation.  
[1, <https://en.wikipedia.org/wiki/Propofol>, 23.12.2019 13:54]
- Propofol is also used in veterinary medicine for anesthesia. [1]
- Synonyms 2,6-bis(1-methylethyl)phenol, 2,6-Diisopropylphenol  
[2, <https://www.drugbank.ca/drugs/DB00818>, 23.12.2019 13:59]
- Product names: Diprivan, Propofol, Propofol Injection, PMS-propofol, Teva-propofol, Propofol Injectable Emulsion, Anesthesia S/I-40 – mixture product, Propoven, Anepol, Anespro Anesvan, Critifol, Disoprivan, Disoprofol / Dormofol, Fresofol, Gobbifol, Hipnolam, Hypro, IV-Pro, Lipuro, Oleo-Lax, Plofed, Profol, Profolen, Propofabb, Propofil, Propogen, Propolipid, Propovan, Propoven, Provive, Rapinovet, Recofol, Safol, Trivam, Troypofol, Unifol [2]

# Challenge

- After changes in production formerly acceptable lipid emulsions with Propofol changed within few months.
- Analytical evaluation by laser diffraction (LD) and dynamic light scattering (DLS) did not reveal differences in the fresh product just after production.
- After several months additional droplet fractions were measured by LD, DLS, but these techniques do not measure the emulsions in the original concentration. Nor do they allow for a prediction.
- The need for an early and reliable product characterization exists to avoid the loss of large product quantities / the recall of batches / the loss of profit.

# Solution



- ✓ Direct characterization of the lipid emulsions with Propofol in original concentration (in contrast to dilute samples for LD, DLS) by accelerated separation.
- ✓ Qualitative and quantitative comparison of acceptable and non-acceptable batches.
- ✓ Complimentary determination of the droplet size distribution for the dilute samples to compare these results with previous analytical approaches.

# Lipid emulsion with Propofol

38 h at gravity

vs.

2 h at higher gravity




38 h real-time separation at gravity, RT, every 10<sup>th</sup> profile shown.

125 minutes accelerated separation at RCA 2300, RT, every 10<sup>th</sup> profile shown.


Identical qualitative separation behaviour allows for quick analysis at high gravity.

# Lipid emulsion with Propofol




125 minutes accelerated separation at RCA 2300, 25 °C, every 10<sup>th</sup> profile shown, 870 nm, photo after 125 min.

# Fast comparison of original samples using creaming velocity, no dilution



50 minutes at RCA 2300, 25 °C, Tracking of creaming front at 10 % transmission

# Droplet size distribution, dilute samples




Droplet size distributions of creaming particles, calculated with customer's material data, are in agreement with the results from other particle sizing methods. Dilution for LUMiSizer 1:20 (m/m).



# Accelerated measurements also allow for prediction of separation at gravity

Creaming velocity = f (RCA) for three emulsions A, B, C



All three emulsions show a linear dependence of their creaming velocity on RCA. The calculation of the creaming velocity at gravity (RCA=1) is possible by linear regression.



# Summary



- ✓ LUMiSizer allows for the fast and direct characterization of the lipid emulsions with Propofol in original concentration by accelerated separation (in contrast to dilute samples for LD, DLS).
- ✓ Higher gravity and earth gravity separation show identical qualitative behaviour.
- ✓ Visual inspection after separation is in agreement with last transmission profiles at 870 nm.
- ✓ Batch release is easy by comparative analysis applying ISO/TR 13097.
- ✓ Complimentary determination of the droplet size distribution for the dilute samples is in agreement with LD/DLS results. Differences between benchmark and acceptable batches are smaller than to failed batches.
- ✓ Comparative and predictive analysis according to ISO/TR 13097 is enabled.