

## Characterizing Gold Nanorods by Utilizing VV and VH Modes in DLS

Zhibin Guo, Hui Ning

Application Research Lab, Bettersize Instruments Ltd.

### Introduction





Figure 1. SEM Images of Gold Nanorods

Some samples exhibit anisotropic shapes such as rods or disks rather than spherical forms. In dynamic light scattering (DLS) measurements, the diffusion behavior of these samples involves not only translational diffusion but also significant rotational diffusion. For these samples, it is essential to perform DLS measurements using both VV and VH modes.

- VV mode uses a vertically polarized laser source and a vertical polarizer in the detection path, collecting only vertically polarized scattered light.
- VH mode uses a vertically polarized laser source and a horizontal polarizer in the detection path, collecting only horizontally polarized scattered light. VH mode is also known as depolarized dynamic light scattering (DDLS).

Combining the results from VV and VH modes allows for the determination of the rotational diffusion coefficient  $D_r$  in anisotropic systems, providing more information about the sample's morphology and motion.

In this application note, the BeNano 180 Zeta Pro is used to measure the translational and rotational diffusion of gold nanorods dispersed in water.



BeNano 180 Zeta Pro

## Experimental

The gold nanorods were measured using the BeNano 180 Zeta Pro with a built-in temperature control system, set to a default test temperature of  $25^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$ .



Figure 2. Multifunctional Micro-volume Cells for VH (left) and VV (right) Modes

Using the micro-volume cells, only  $16\ \mu\text{L}$  of the sample is needed, significantly saving valuable sample. The samples were measured at least three times to ensure reproducibility and to calculate the standard deviation of the results.

## Results and Discussion

VV and VH mode tests were performed at a  $90^{\circ}$  angle. The DLS results are shown below:




Figure 3.  $90^{\circ}$  VV Mode Correlation Functions of Multiple Tests




Figure 4.  $90^{\circ}$  VH Mode Correlation Functions of Multiple Tests

The derived formula from literature shows that the light intensity correlation function for rod-shaped samples under VV and VH test conditions is:

$$g_2^{VV} - 1 = \beta[A^2 \exp\{-\Gamma t\} + 2AB \exp\{-(\Gamma + \Delta/2)t\} + B^2 \exp\{-(\Gamma + \Delta)t\}]$$

$$g_2^{VH} - 1 = \beta' \exp\{-(\Gamma + \Delta)t\}$$

where  $\Gamma = 2D_T q^2$ ,  $\Delta = 12D_r$ ,  $q$  is the scattering vector. In the VV mode, the correlation function includes contributions from both translational and rotational diffusion, leading to multiple exponential decays. In the VH mode, the correlation function shows a single exponential decay. This matches the observed correlation functions in Figures 3 and 4.

Based on theoretical derivations, the field intensity correlation functions for VV and VH modes are as follows:

$$g_1^{VV} = A \exp\{-D_T q^2 t\} + B \exp\{-(D_T q^2 + 6D_r)t\}$$

$$g_1^{VH} = \exp\{-(D_T q^2 + 6D_r)t\}$$




Figure 5. Double Exponential Fitting of  $g_1^{VV}$  Correlation Function

By calculating the natural logarithm of  $g_1^{VH}$  and using polynomial fitting, one can obtain  $D_T q^2 + 6D_r$ . Double exponential fitting of  $g_1^{VV}$  yields both  $D_T q^2 + 6D_r$  and  $D_T q^2$ . The results are summarized in the table below:

Table 1. Fitting results from VV and VH modes

Item	Method	Result
$D_T q^2 + 6D_r$ (1/s)	Polynomial fit of $g_1^{VH}$	66935.67
$D_T q^2 + 6D_r$ (1/s)	Fast decay in VV mode	63260.50
$D_T$ ( $\text{nm}^2/\text{ms}$ )	Slow decay in VV mode	2941.07
$D_r$ (1/ms)	Combined fast and slow decay in VV mode	10.58

Table 1 shows that the results from VV and VH modes are consistent, ultimately providing the rotational diffusion coefficient  $D_r$  for the gold nanorods.

## Conclusion

The combined use of VV and VH modes in dynamic light scattering (DLS) tests allows for comprehensive analysis of anisotropic particles like gold nanorods. The BeNano 180 Zeta Pro effectively measures both translational and rotational diffusion, providing detailed insights into the particles' morphology and motion. This approach ensures accurate characterization and enhances understanding of such complex systems.

# Bettersize

BETTER PARTICLE SIZE SOLUTIONS

[www.bettersizeinstruments.com](http://www.bettersizeinstruments.com)

[info@bettersize.com](mailto:info@bettersize.com)

## **Bettersize Instruments Ltd.**

**Address:** No. 9, Ganquan Road, Jinquan Industrial Park,

Dandong, Liaoning, China

**Postcode:** 118009

**Tel:** +86-755-26926582

**Fax:** +86-415-6170645

## **Bettersize Inc.**

**Address:** Suite K-2, 3188 Airway Ave, Costa Mesa, CA 92626,

United States

**Tel:** +1 833-699-7493 (SIZE)

Download Our Application Notes:



Visit Our BeNano 180 Zeta Pro Site:



**Disclaimer:** By using or accessing the brochure, you agree with the Disclaimer without any qualification or limitation. Diligent care has been used to ensure that the information in this brochure is accurate, Bettersize Instruments Ltd. shall not be liable for errors contained herein or for damages in connection with the use of this material. The information on this brochure is presented as general information and no representation or warranty is expressly or impliedly given as to its accuracy, completeness or correctness. It does not constitute part of a legal offer or contract. Bettersize Instruments Ltd. reserves the right to modify, alter, add and delete the content outlined in the brochure without prior notice and without any subsequent liability to the company.

Copyright: © 2024 Bettersize Instruments Ltd. | All Rights Reserved