

METHOD 3550B

ULTRASONIC EXTRACTION

See Disclaimer. See manufacturer's specifications for operational settings.

1.0 SCOPE AND APPLICATION

1.1 Method 3550 is a procedure for extracting nonvolatile and semivolatile organic compounds from solids such as soils, sludges, and wastes. The ultrasonic process ensures intimate contact of the sample matrix with the extraction solvent.

1.2 The method is divided into two sections, based on the expected concentration of organics in the sample. The low concentration method (individual organic components of less than or equal to 20 mg/kg) uses a larger sample size and a more rigorous extraction procedure (lower concentrations are more difficult to extract). The medium/high concentration method (individual organic components of greater than 20 mg/kg) is much simpler and therefore faster.

1.3 It is highly recommended that the extracts be cleaned up prior to analysis. See Chapter Four (Cleanup), Sec. 4.2.2, for applicable methods.

1.4 Ultrasonic extraction is not as rigorous as other extraction methods for soils/solids. Therefore, it is critical that the method (including the manufacturer's instructions) be followed explicitly in order to achieve the maximum extraction efficiency. See Sec. 7.0 for the critical aspects of the extraction procedure.

1.5 EPA has not validated Method 3550 for the extraction of organophosphorous compounds from solid matrices. In addition, there are concerns that the ultrasonic energy may lead to breakdown of some organophosphorous compounds (see Reference 3). As a result, this extraction technique should not be used for organophosphorous compounds without extensive validation on real-world samples. Such studies should assess the precision, accuracy, ruggedness, and sensitivity of the technique relative to the appropriate regulatory limits or project-specific concentrations of interest.

1.6 This method is restricted to use by or under the supervision of trained analysts. Each analyst must demonstrate the ability to generate acceptable results with this method.

1.7 This method is not appropriate for applications where high extraction efficiencies of analytes at very low concentrations is necessary (e.g., demonstration of effectiveness of corrective action).

2.0 SUMMARY OF METHOD

2.1 Low concentration method - A 30-g sample is mixed with anhydrous sodium sulfate to form a free-flowing powder. This is solvent extracted three times using ultrasonic extraction. The extract is separated from the sample by vacuum filtration or centrifugation. The extract is ready for cleanup and/or analysis following concentration.

2.2 Medium/high concentration method - A 2-g sample is mixed with anhydrous sodium sulfate to form a free-flowing powder. This is solvent extracted once using ultrasonic extraction. A portion of the extract is removed for cleanup and/or analysis.

3.0 INTERFERENCES

Refer to Method 3500.

4.0 APPARATUS AND MATERIALS

4.1 Apparatus for grinding dry waste samples.

4.2 Ultrasonic preparation - A horn-type device equipped with a titanium tip, or a device that will give equivalent performance, shall be used.

4.2.1 Ultrasonic Disrupter - The disrupter must have a minimum power wattage of 300 watts, with pulsing capability. A device designed to reduce the cavitation sound is recommended. Follow the manufacturers instructions for preparing the disrupter for extraction of samples with low and medium/high concentration.

4.2.2 Use a 3/4" horn for the low concentration method and a 1/8" tapered microtip attached to a 1/2" horn for the medium/high concentration method.

4.3 Sonobox - Recommended with above disrupters for decreasing cavitation sound (Heat Systems - Ultrasonics, Inc., Model 432B or equivalent).

4.4 Apparatus for determining percent dry weight.

4.4.1 Drying oven - capable of maintaining 105°C.

4.4.2 Desiccator.

4.4.3 Crucibles - Porcelain or disposable aluminum.

4.5 Pasteur glass pipets - 1-mL, disposable.

4.6 Beakers - 400-mL.

4.7 Vacuum or pressure filtration apparatus.

4.7.1 Buchner funnel.

4.7.2 Filter paper - Whatman No. 41 or equivalent.

4.8 Kuderna-Danish (K-D) apparatus.

4.8.1 Concentrator tube - 10-mL, graduated (Kontes K-570050-1025 or equivalent). A ground-glass stopper is used to prevent evaporation of extracts.

4.8.2 Evaporation flask - 500-mL (Kontes K-570001-500 or equivalent). Attach to concentrator tube with springs, clamps, or equivalent.

4.8.3 Snyder column - Three-ball macro (Kontes K-503000-0121 or equivalent).

4.8.4 Snyder column - Two-ball micro (Kontes K-569001-0219 or equivalent).

4.8.5 Springs - 1/2 inch (Kontes K-662750 or equivalent).

NOTE: The following glassware is recommended for the purpose of solvent recovery during the concentration procedures requiring the use of Kuderna-Danish evaporative concentrators. Incorporation of this apparatus may be required by State or local municipality regulations that govern air emissions of volatile organics. EPA recommends the incorporation of this type of reclamation system as a method to implement an emissions reduction program. Solvent recovery is a means to conform with waste minimization and pollution prevention initiatives.

4.9 Solvent vapor recovery system (Kontes K-545000-1006 or K-547300-0000, Ace Glass 6614-30, or equivalent).

4.10 Boiling chips - Solvent-extracted, approximately 10/40 mesh (silicon carbide or equivalent).

4.11 Water bath - Heated, with concentric ring cover, capable of temperature control ($\pm 5^{\circ}\text{C}$). The bath should be used in a hood.

4.12 Balance - Top-loading, capable of accurately weighing to the nearest 0.01 g.

4.13 Vials - 2-mL, for GC autosampler, with polytetrafluoroethylene (PTFE)-lined screw caps or crimp tops.

4.14 Glass scintillation vials - 20-mL, with PTFE-lined screw caps.

4.15 Spatula - Stainless steel or PTFE.

4.16 Drying column - 20-mm ID Pyrex chromatographic column with Pyrex glass wool at bottom.

NOTE: Fritted glass discs are difficult to decontaminate after highly contaminated extracts have been passed through. Columns without frits may be purchased. Use a small pad of Pyrex glass wool to retain the adsorbent. Prewash the glass wool pad with 50 mL of acetone followed by 50 mL of elution solvent prior to packing the column with adsorbent.

4.17 Syringe - 5-mL.

5.0 REAGENTS

5.1 Reagent grade inorganic chemicals shall be used in all tests. Unless otherwise specified, it is intended that all inorganic reagents shall conform to the specifications of the Committee on

Analytical Reagents of the American Chemical Society, where such specifications are available. Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination.

5.2 Organic-free reagent water. All references to water in this method refer to organic-free reagent water, as defined in Chapter One.

5.3 Sodium sulfate (granular, anhydrous), Na_2SO_4 . Purify by heating at 400°C for 4 hours in a shallow tray, or by precleaning the sodium sulfate with methylene chloride. If the sodium sulfate is precleaned with methylene chloride, a method blank must be analyzed, demonstrating that there is no interference from the sodium sulfate.

5.4 Extraction solvents

Samples should be extracted using a solvent system that gives optimum, reproducible recovery of the analytes of interest from the sample matrix. Table 1 provides recovery data for selected semivolatile organic compounds extracted from an NIST SRM. The following sections provide guidance on the choice of solvents for various classes of analytes. All solvents must be pesticide quality or equivalent.

5.4.1 Semivolatile organics may be extracted with acetone/methylene chloride (1:1, v/v), $\text{CH}_3\text{COCH}_3/\text{CH}_2\text{Cl}_2$ or acetone/hexane (1:1, v/v), $\text{CH}_3\text{COCH}_3/\text{C}_6\text{H}_{14}$.

5.4.2 Organochlorine pesticides may be extracted with acetone/hexane (1:1, v/v), $\text{CH}_3\text{COCH}_3/\text{C}_6\text{H}_{14}$ or acetone/methylene chloride (1:1, v/v), $\text{CH}_3\text{COCH}_3/\text{CH}_2\text{Cl}_2$.

5.4.3 PCBs may be extracted with acetone/hexane (1:1, v/v), $\text{CH}_3\text{COCH}_3/\text{C}_6\text{H}_{14}$, acetone/methylene chloride (1:1, v/v), $\text{CH}_3\text{COCH}_3/\text{CH}_2\text{Cl}_2$ or hexane, C_6H_{14} .

5.4.4 Other solvent systems may be employed, provided that the analyst can demonstrate adequate performance for the analytes of interest in the sample matrix (see Method 3500, Sec. 8.0).

5.5 Exchange solvents - All solvents must be pesticide quality or equivalent.

5.5.1 Hexane, C_6H_{14} .

5.5.2 2-Propanol, $(\text{CH}_3)_2\text{CHOH}$.

5.5.3 Cyclohexane, C_6H_{12} .

5.5.4 Acetonitrile, CH_3CN .

5.5.5 Methanol, CH_3OH .

6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

See the introductory material to this chapter, Organic Analytes, Sec. 4.1.

7.0 PROCEDURE

As noted in Sec. 1.4, ultrasonic extraction is not as rigorous a method as other extraction methods for soils/solids. Therefore, it is critical that the method be followed explicitly (including the manufacturer's instructions) to achieve the maximum extraction efficiency. At a minimum, successful use of this technique requires that:

- The extraction device must have a minimum of 300 watts of power and be equipped with appropriate size disrupter horns (see Sec. 4.2).
- The horn must be properly maintained, including tuning according to the manufacturer's instructions prior to use, and inspection of the horn tip for excessive wear.
- The samples must be properly prepared by thorough mixing with sodium sulfate so that it forms a free-flowing powder prior to the addition of the solvent.
- The extraction horns used for the low concentration and high concentration protocols (Sec. 7.3 and Sec. 7.4, respectively) are not interchangeable. Results indicate that the use of the 3/4" horn is inappropriate for the high concentration method, particularly for extraction of very non-polar organic compounds such as PCBs, which are strongly adsorbed to the soil matrix.
- Three extractions are performed with the appropriate solvent, the extraction is performed in the specified pulse mode, and the horn tip is positioned just below the surface of the solvent yet above the sample.
- Very active mixing of the sample and the solvent must occur when the ultrasonic pulse is activated. The analyst must observe such mixing at some point during the extraction process.

7.1 Sample handling

7.1.1 Sediment/soil samples - Decant and discard any water layer on a sediment sample. Mix sample thoroughly, especially composited samples. Discard any foreign objects such as sticks, leaves, and rocks.

7.1.2 Waste samples - Samples consisting of multiple phases must be prepared by the phase separation method in Chapter Two before extraction. This extraction procedure is for solids only.

7.1.3 Dry waste samples amenable to grinding - Grind or otherwise subdivide the waste so that it either passes through a 1-mm sieve or can be extruded through a 1-mm hole. Introduce sufficient sample into the grinding apparatus to yield at least 10 g after grinding.

7.1.4 Gummy, fibrous, or oily materials not amenable to grinding should be cut, shredded, or otherwise reduced in size to allow mixing and maximum exposure of the sample surfaces for the extraction. The addition of anhydrous sodium sulfate to the sample (1:1) may make the mixture amenable to grinding.

7.2 Determination of percent dry weight - When sample results are to be calculated on a dry weight basis, a second portion of sample should be weighed out at the same time as the portion used for analytical determination.

WARNING: The drying oven should be contained in a hood or vented. Significant laboratory contamination may result from drying a heavily contaminated sample.

Immediately after weighing the sample for extraction, weigh 5-10 g of the sample into a tared crucible. Dry this aliquot overnight at 105°C. Allow to cool in a desiccator before weighing. Calculate the % dry weight as follows:

$$\% \text{ dry weight} = \frac{\text{g of dry sample}}{\text{g of sample}} \times 100$$

7.3 Extraction method for samples expected to contain low concentrations of organics and pesticides (less than or equal to 20 mg/kg):

7.3.1 The following steps should be performed rapidly to avoid loss of the more volatile extractables.

7.3.1.1 Weigh approximately 30 g of sample into a 400-mL beaker. Record the weight to the nearest 0.1 g.

7.3.1.2 Nonporous or wet samples (gummy or clay type) that do not have a free-flowing sandy texture must be mixed with 60 g of anhydrous sodium sulfate, using a spatula. If required, more sodium sulfate may be added. After addition of sodium sulfate, the sample should be free flowing.

7.3.1.3 Add 1.0 mL of the surrogate standard solution to all samples, spiked samples, QC samples, and blanks. Consult Method 3500, Secs. 5.0 and 8.0 for guidance on the appropriate choice of surrogate compounds and concentrations.

7.3.1.4 For the sample in each batch selected for spiking, add 1.0 mL of the matrix spiking solution. Consult Method 3500, Secs. 5.0 and 8.0 for guidance on the appropriate choice of matrix spiking compounds and concentrations.

7.3.1.5 If gel permeation cleanup (Method 3640) is to be employed, the analyst should either add twice the volume of the surrogate spiking solution (and matrix spiking solution, where applicable), or concentrate the final extract to half the normal volume, to compensate for the half of the extract that is lost due to loading of the GPC column.

7.3.1.6 Immediately add 100 mL of the appropriate/recommended extraction solvent or solvent mixture (see Sec. 5.4 and Table 1).

7.3.2 Place the bottom surface of the tip of the #207 (or equivalent) 3/4 inch disrupter horn about 1/2 inch below the surface of the solvent, but above the sediment layer.

NOTE: Be sure the horn is properly tuned according to the manufacturer's instructions.

7.3.3 Extract ultrasonically for 3 minutes, with output control knob set at 10 (full power) and with mode switch on Pulse (pulsing energy rather than continuous energy) and percent-duty cycle knob set at 50% (energy on 50% of time and off 50% of time). Do not use microtip probe.

7.3.4 Decant the extract and filter it through Whatman No. 41 filter paper (or equivalent) in a Buchner funnel that is attached to a clean 500-mL filtration flask. Alternatively, decant the extract into a centrifuge bottle and centrifuge at low speed to remove particles.

7.3.5 Repeat the extraction two or more times with two additional 100 mL portions of solvent. Decant off the solvent after each ultrasonic extraction. On the final ultrasonic extraction, pour the entire sample into the Buchner funnel and rinse with extraction solvent. Apply a vacuum to the filtration flask, and collect the solvent extract. Continue filtration until all visible solvent is removed from the funnel, but do not attempt to completely dry the sample, as the continued application of a vacuum may result in the loss of some analytes. Alternatively, if centrifugation is used in Sec. 7.3.4, transfer the entire sample to the centrifuge bottle. Centrifuge at low speed, and then decant the solvent from the bottle.

7.3.6 Assemble a Kuderna-Danish (K-D) concentrator (if necessary) by attaching a 10-mL concentrator tube to a 500-mL evaporator flask. Attach the solvent vapor recovery glassware (condenser and collection device) to the Snyder column of the Kuderna-Danish apparatus following manufacturer's instructions. Transfer filtered extract to a 500-mL evaporator flask and proceed to the next section.

7.3.7 Add one to two clean boiling chips to the evaporation flask, and attach a three-ball Snyder column. Prewet the Snyder column by adding about 1 mL methylene chloride to the top. Place the K-D apparatus on a hot water bath (80 - 90°C) so that the concentrator tube is partially immersed in the hot water and the entire lower rounded surface of the flask is bathed with hot vapor. Adjust the vertical position of the apparatus and the water temperature, as required, to complete the concentration in 10 - 15 min. At the proper rate of distillation the balls of the column will actively chatter, but the chambers will not flood with condensed solvent. When the apparent volume of liquid reaches 1 mL, remove the K-D apparatus and allow it to drain and cool for at least 10 min.

7.3.8 If a solvent exchange is required (as indicated in Table 1), momentarily remove the Snyder column, add 50 mL of the exchange solvent and a new boiling chip, and re-attach the Snyder column. Concentrate the extract as described in Sec. 7.3.10, raising the temperature of the water bath, if necessary, to maintain proper distillation. When the apparent volume again reaches 1 - 2 mL, remove the K-D apparatus and allow it to drain and cool for at least 10 minutes.

7.3.9 Remove the Snyder column and rinse the flask and its lower joints into the concentrator tube with 1 - 2 mL of methylene chloride or exchange solvent. If sulfur crystals are a problem, proceed to Method 3660 for cleanup. The extract may be further concentrated by using the technique outlined in Sec. 7.3.10 or adjusted to 10.0 mL with the solvent last used.

7.3.10 If further concentration is indicated in Table 1, either micro Snyder column technique (Sec. 7.3.10.1) or nitrogen blowdown technique (Sec. 7.3.10.2) may be used to adjust the extract to the final volume required.

7.3.10.1 Micro Snyder column technique

7.3.10.1.1 Add a clean boiling chip and attach a two-ball micro Snyder column to the concentrator tube. Prewet the column by adding approximately 0.5 mL of methylene chloride or exchange solvent through the top. Place the apparatus in the hot water bath. Adjust the vertical position and the

water temperature, as required, to complete the concentration in 5-10 minutes. At the proper rate of distillation the balls of the column will actively chatter, but the chambers will not flood. When the liquid reaches an apparent volume of approximately 0.5 mL, remove the apparatus from the water bath and allow to drain and cool for at least 10 minutes. Remove the micro Snyder column and rinse its lower joint with approximately 0.2 mL of appropriate solvent and add to the concentrator tube. Adjust the final volume to the volume required for cleanup or for the determinative method (see Table 1).

7.3.10.2 Nitrogen blowdown technique

7.3.10.2.1 Place the concentrator tube in a warm water bath (approximately 35°C) and evaporate the solvent volume to the required level using a gentle stream of clean, dry nitrogen (filtered through a column of activated carbon).

CAUTION: Do not use plasticized tubing between the carbon trap and the sample, since it may introduce contaminants.

7.3.10.2.2 The internal wall of the tube must be rinsed down several times with the appropriate solvent during the operation. During evaporation, the solvent level in the tube must be positioned to prevent water from condensing into the sample (i.e., the solvent level should be below the level of the water bath). Under normal operating conditions, the extract should not be allowed to become dry.

CAUTION: When the volume of solvent is reduced below 1 mL, semivolatile analytes may be lost.

7.4 Extraction method for samples expected to contain high concentrations of organics (greater than 20 mg/kg):

7.4.1 Transfer approximately 2 g (record weight to the nearest 0.1 g) of sample to a 20-mL vial. Wipe the mouth of the vial with a tissue to remove any sample material. Record the exact weight of sample taken. Cap the vial before proceeding with the next sample to avoid any cross contamination.

7.4.2 Add 2 g of anhydrous sodium sulfate to sample in the 20-mL vial and mix well.

7.4.3 Surrogates are added to all samples, spikes, and blanks (see Method 3500 for details on the surrogate spiking solution and on the matrix spike solution).

7.4.3.1 Add 1.0 mL of surrogate spiking solution to sample mixture.

7.4.3.2 For the sample in each analytical batch selected for spiking, add 1.0 mL of the matrix spiking standard.

7.4.3.3 If gel permeation cleanup (Method 3640) is to be employed, the analyst should either add twice the volume of the surrogate spiking solution (and matrix spiking solution, where applicable), or concentrate the final extract to half the normal volume, to compensate for the half of the extract that is lost due to loading of the GPC column.

7.4.4 Immediately add whatever volume of solvent is necessary to bring the final volume to 10.0 mL considering the added volume of surrogates and matrix spikes. Disrupt the sample with the 1/8 in. tapered microtip ultrasonic probe for 2 minutes at output control setting 5 and with mode switch on pulse and percent duty cycle at 50%. Extraction solvents are:

7.4.4.1 For nonpolar compounds (i.e., organochlorine pesticides and PCBs), use hexane or appropriate solvent.

7.4.4.2 For other semivolatile organics, use methylene chloride.

7.4.5 Loosely pack a disposable Pasteur pipette with 2 to 3 cm of glass wool. Filter the sample extract through the glass wool and collect the extract in a suitable container. The entire 10 mL of extraction solvent cannot be recovered from the sample. Therefore, the analyst should collect a volume appropriate for the sensitivity of the determinative method. For instance, for methods that do not require that the extract be concentrated further (e.g., Method 8081 typically employs a final extract volume of 10 mL), the extract may be collected in a scintillation vial or other sealable container. For extracts that will require further concentration, it is advisable to collect a standard volume for all such samples in order to simplify the calculation of the final sample results. For instance, collect 5.0 mL of extract in a clean concentrator tube. This volume represents exactly half of the total volume of the original sample extract. As necessary, account for the "loss" of half of the extract in the final sample calculations, or concentrate the final extract to one-half the nominal final volume (e.g., 0.5 mL vs. 1.0 mL) to compensate for the loss.

7.4.6 The extract is ready for cleanup or analysis, depending on the extent of interfering co-extractives.

7.5 If analysis of the extract will not be performed immediately, stopper the concentrator tube and refrigerate. If the extract will be stored longer than 2 days, it should be transferred to a vial with a PTFE-lined cap and labeled appropriately.

8.0 QUALITY CONTROL

8.1 Any reagent blanks, matrix spike, and replicate samples should be subjected to exactly the same analytical procedures as those used on actual samples.

8.2 Refer to Chapter One for specific quality control procedures and Method 3500 for extraction and sample preparation procedures.

9.0 METHOD PERFORMANCE

Refer to the determinative methods for performance data.

10.0 REFERENCES

1. U.S. EPA, Interlaboratory Comparison Study: Methods for Volatile and Semi-Volatile Compounds, Environmental Monitoring Systems Laboratory, Office of Research and Development, Las Vegas, NV, EPA 600/4-84-027, 1984.

2. Christopher S. Hein, Paul J. Marsden, Arthur S. Shurtleff, "Evaluation of Methods 3540 (Soxhlet) and 3550 (Sonication) for Evaluation of Appendix IX Analytes from Solid Samples", S-CUBED, Report for EPA Contract 68-03-33-75, Work Assignment No. 03, Document No. SSS-R-88-9436, October 1988.
3. Kotronarou, A., et al., "Decomposition of Parathion in Aqueous Solution by Ultrasonic Irradiation," *ES&T*, 1992, Vol. 26, 1460-1462.

TABLE 1
EFFICIENCIES OF VARIOUS EXTRACTION SOLVENT SYSTEMS^a

| Compound | CAS No. ^b | ABN ^c | Solvent System ^d | | | | | | | | | |
|-------------------------------|----------------------|------------------|-----------------------------|------|------|-----|-------|-----|------|------|------|-----|
| | | | A | | B | | C | | D | | E | |
| | | | %R | SD | %R | SD | %R | SD | %R | SD | %R | SD |
| 4-Bromophenyl phenyl ether | 101-55-3 | N | 64.2 | 6.5 | 56.4 | 0.5 | 86.7 | 1.9 | 84.5 | 0.4 | 73.4 | 1.0 |
| 4-Chloro-3-methylphenol | 59-50-7 | A | 66.7 | 6.4 | 74.3 | 2.8 | 97.4 | 3.4 | 89.4 | 3.8 | 84.1 | 1.6 |
| Bis(2-chloroethoxy)methane | 111-91-1 | N | 71.2 | 4.5 | 58.3 | 5.4 | 69.3 | 2.4 | 74.8 | 4.3 | 37.5 | 5.8 |
| Bis(2-chloroethyl) ether | 111-44-4 | N | 42.0 | 4.8 | 17.2 | 3.1 | 41.2 | 8.4 | 61.3 | 11.7 | 4.8 | 1.0 |
| 2-Chloronaphthalene | 91-58-7 | N | 86.4 | 8.8 | 78.9 | 3.2 | 100.8 | 3.2 | 83.0 | 4.6 | 57.0 | 2.2 |
| 4-Chlorophenyl phenyl ether | 7005-72-3 | N | 68.2 | 8.1 | 63.0 | 2.5 | 96.6 | 2.5 | 80.7 | 1.0 | 67.8 | 1.0 |
| 1,2-Dichlorobenzene | 95-50-1 | N | 33.3 | 4.5 | 15.8 | 2.0 | 27.8 | 6.5 | 53.2 | 10.1 | 2.0 | 1.2 |
| 1,3-Dichlorobenzene | 541-73-1 | N | 29.3 | 4.8 | 12.7 | 1.7 | 20.5 | 6.2 | 46.8 | 10.5 | 0.6 | 0.6 |
| Diethyl phthalate | 84-66-2 | N | 24.8 | 1.6 | 23.3 | 0.3 | 121.1 | 3.3 | 99.0 | 4.5 | 94.8 | 2.9 |
| 4,6-Dinitro- <i>o</i> -cresol | 534-52-1 | A | 66.1 | 8.0 | 63.8 | 2.5 | 74.2 | 3.5 | 55.2 | 5.6 | 63.4 | 2.0 |
| 2,4-Dinitrotoluene | 121-14-2 | N | 68.9 | 1.6 | 65.6 | 4.9 | 85.6 | 1.7 | 68.4 | 3.0 | 64.9 | 2.3 |
| 2,6-Dinitrotoluene | 606-20-2 | N | 70.0 | 7.6 | 68.3 | 0.7 | 88.3 | 4.0 | 65.2 | 2.0 | 59.8 | 0.8 |
| Heptachlor epoxide | 1024-57-3 | N | 65.5 | 7.8 | 58.7 | 1.0 | 86.7 | 1.0 | 84.8 | 2.5 | 77.0 | 0.7 |
| Hexachlorobenzene | 118-74-1 | N | 62.1 | 8.8 | 56.5 | 1.2 | 95.8 | 2.5 | 89.3 | 1.2 | 78.1 | 4.4 |
| Hexachlorobutadiene | 87-68-3 | N | 55.8 | 8.3 | 41.0 | 2.7 | 63.4 | 4.1 | 76.9 | 8.4 | 12.5 | 4.6 |
| Hexachlorocyclopentadiene | 77-47-4 | N | 26.8 | 3.3 | 19.3 | 1.8 | 35.5 | 6.5 | 46.6 | 4.7 | 9.2 | 1.7 |
| Hexachloroethane | 67-72-1 | N | 28.4 | 3.8 | 15.5 | 1.6 | 31.1 | 7.4 | 57.9 | 10.4 | 1.4 | 1.2 |
| 5-Nitro- <i>o</i> -toluidine | 99-55-8 | B | 52.6 | 26.7 | 64.6 | 4.7 | 74.7 | 4.7 | 27.9 | 4.0 | 34.0 | 4.0 |
| Nitrobenzene | 98-95-3 | N | 59.8 | 7.0 | 38.7 | 5.5 | 46.9 | 6.3 | 60.6 | 6.3 | 13.6 | 3.2 |
| Phenol | 108-95-2 | A | 51.6 | 2.4 | 52.0 | 3.3 | 65.6 | 3.4 | 65.5 | 2.1 | 50.0 | 8.1 |
| 1,2,4-Trichlorobenzene | 120-82-1 | N | 66.7 | 5.5 | 49.9 | 4.0 | 73.4 | 3.6 | 84.0 | 7.0 | 20.0 | 3.2 |

^a Percent recovery of analytes spiked at 200 mg/kg into NIST sediment SRM 1645

^b Chemical Abstracts Service Registry Number

^c Compound Type: A = Acid, B = Base, N = neutral

^d Solvent system A = Methylene chloride
 Solvent system B = Methylene chloride/Acetone (1/1)
 Solvent system C = Hexane/Acetone (1/1)
 Solvent system D = Methyl t-butyl ether
 Solvent system E = Methyl t-butyl ether/Methanol (2/1)

TABLE 2
SPECIFIC EXTRACTION CONDITIONS FOR VARIOUS DETERMINATIVE METHODS


| Determinative method | Extraction pH | Exchange solvent required for analysis | Exchange solvent required for cleanup | Volume of extract required for cleanup (mL) | Final extract volume for analysis (mL) ^a |
|----------------------|---------------|--|---------------------------------------|---|---|
| 8041 | as received | 2-propanol | hexane | 1.0 | 1.0, 0.5 ^b |
| 8061 | as received | hexane | hexane | 2.0 | 10.0 |
| 8070 | as received | methanol | methylene chloride | 2.0 | 10.0 |
| 8081 | as received | hexane | hexane | 10.0 | 10.0 |
| 8091 | as received | hexane | hexane | 2.0 | 1.0 |
| 8100 | as received | none | cyclohexane | 2.0 | 1.0 |
| 8111 | as received | hexane | hexane | 2.0 | 10.0 |
| 8121 | as received | hexane | hexane | 2.0 | 1.0 |
| 8141 | as received | hexane | hexane | 10.0 | 10.0 |
| 8270 ^c | as received | none | - | - | 1.0 |
| 8310 | as received | acetonitrile | - | - | 1.0 |
| 8321 | as received | methanol | - | - | 1.0 |
| 8325 | as received | methanol | - | - | 1.0 |
| 8410 | as received | methylene chloride | methylene chloride | 10.0 | 0.0 (dry) |

^a For methods where the suggested final extract volume is 10.0 mL, the volume may be reduced to as low as 1.0 mL to achieve lower detection limits.

^b Phenols may be analyzed by Method 8041, using a 1.0-mL 2-propanol extract by GC/FID. Method 8041 also contains an optional derivatization procedure for phenols which results in a 0.5-mL hexane extract to be analyzed by GC/ECD.

^c The specificity of GC/MS may make cleanup of the extracts unnecessary. Refer to Method 3600 for guidance on the cleanup procedures available if required.

METHOD 3550B
ULTRASONIC EXTRACTION



METHOD 3550B
continued

