



# xcell

## Optimized virus production HEK ViP media family

- For transfection and production of proteins and viral vectors (e.g. AAV and lentivirus) in HEK and other human cell lines
- Comprehensive solution:** thawing, expanding, transfecting and producing can all be performed in the same medium, regardless of the scale
- HEK ViP NB** and **HEK ViP NX:** basic and high nutrient levels to fit each specific application and cell line
- Chemically defined, protein-free, AOF, w/o growth factor



### Xcell AG

Alte Verler Strasse 1  
33689 Bielefeld  
Germany

+49 (0)521 96989-200

+49 (0)521 96989-201

info@xcell.de

www.xcell.de